

Njord Micro-Organism Test – Liqui-Cult

Liqui-Cult™ test kits are used to detect microbial growth in fluids, in the laboratory and in the workplace. Areas of application included metalworking fluids, circulating cooling and process water, hydrocarbons, and other liquids such as cleaners, cosmetics, liquid soaps and detergents, etc.


Liqui-Cult™ testing determines the type and degree of infestation usually measured by the more expensive plate count system. To use, simply fill the 5cc syringe (included in your kit) with the fluid to be tested. Inject the complete 5cc s through the rubber cap of test bottle containing the liquid agar and shake it vigorously. Incubate for 30 hrs. for bacterial growth or 72 hrs. for fungal growth and compare it to the appropriate photo charts to determine the degree of contamination.

DIRECTIONS FOR USING THE **Liqui-Cult™** TEST KIT

- 1** Remove center seal from metal cap of test bottle.
- 2** Assemble syringe by first removing plunger nozzle guard. Take cap, not guard, off needle and press into plunger nozzle, leaving needle guard in place until ready to use.
- 3** Remove needle guard and insert tip into liquid to be tested. Draw plunger back, filling syringe with complete 5cc s of fluid.
- 4** Inject needle thru rubber cap on bottle and push down on syringe plunger, emptying contents completely into it. *Remove syringe from bottle, replacing needle guard and discard according to local, state and federal guidelines.* Shake bottle vigorously to mix fluids.
- 5** Incubate bottle at room temperature for 30 hours for microbial growth or 72 hours for fungal growth.
- 6** To measure: Bacterial growth...after 30 hours, compare same bottle to color density chart on reverse side of these instructions.
Fungal growth...after a full 72 hours, compare the test bottle to the density levels indicated on that chart.
- 7** If it is necessary to dilute samples to be tested, for reasons of viscosity, etc., the following procedure is recommended:
A sample may be diluted by adding 1ml to 99ml of boiled, cooled tap water in a clean container. Add 1ml of the test sample to the diluent, cap the container and shake vigorously. Then proceed with the above directions. Take the dilution factor into account when estimating growth. For example, a 10⁶ count from a 1/100 dilution of fluid would indicate 10⁷ organisms per ml in the original sample.
- 8** **STORAGE:** **Liqui-Cult™** should be stored at room temperature. (68-75° degrees F) Protect from light. Need not be refrigerated.
- 9** **DISPOSAL OF USED MATERIALS:** Under a fume hood fill the bottle with bleach...let stand for 24 hours. Then discard in an appropriate manner. *Needles and syringes should be disposed of within local, state and federal guidelines.*


BACTERIAL GROWTH:

This chart is to be used for comparison 30 hours after sample is injected.



FUNGAL GROWTH:

This chart is to be used for comparison 72 hours after sample is injected.



If you have any questions, please contact your nearest Njord Filtration distributor or Njord Filtration directly.